## Nucleis Mag Extractor LIBERTY

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## Nucleis Mag Extractor LIBERTY

### Shipping and storage conditions

Nucleis Mag Extractor LIBERTY is shipped at room temperature. On arrival, the kit should be stored at +4°C. These components remain stable until the indicated expiry date.

### Composition

This kit, for complex matrices such as blood, cerebrospinal fluid or semen, enables up to 14 samples to be extracted from an initial sample volume of 200 µL.

The different products contained in the kit are described in the table below:

| Components        | Volume 14 reactions kit |
|-------------------|-------------------------|
| PK (Proteinase K) | 0.2 mL                  |
| PVG + RNACarrier  | 2.4 mL                  |
| SP                | 0.8 mL                  |
| BMB               | 2 mL                    |
| WASH 1            | 3.5 mL                  |
| WASH 2            | 3.5 mL                  |
| WASH 3            | 3.5 mL                  |
| Elution buffer    | 2 mL                    |

### Equipment required

For manual use of this kit, the equipment (not supplied in the kit) listed below is required:

- Adjustable 10 µL and 200 µL micropipettes
- Sterile low retention filter cones for micropipettes
- Deepwell 96 wells
- Adhesive aluminium foil
- Incubator up to 56°C
- Magnetic support for deepwell 96 wells

## Manual protocol

Depending on the matrix, it is sometimes necessary to carry out a pre-treatment before starting the extraction. The file listing the pre-treatments for each matrix can be found on our website.

1. In a deepwell plate, insert 200  $\mu$ L of sample, add 170  $\mu$ L of PVG + RNACarrier and 10  $\mu$ L of PK. Mix by pipetting gently.
2. Cover the plate with adhesive aluminium foil and incubate for 45 minutes at 56°C.
3. Remove the aluminium foil. Add 40  $\mu$ L of SP and 125  $\mu$ L of BMB, mix by pipetting gently.
4. Incubate for 5 minutes at room temperature.
5. Place the deepwell plate on the magnetic support and remove the supernatant without touching the magnetic beads.
6. Remove the plate from the support, add 200  $\mu$ L of WASH 1 and resuspend the beads.
7. Place the deepwell plate on the magnetic support and remove the supernatant without touching the magnetic beads.
8. Remove the plate from the support, add 200  $\mu$ L of WASH 2 and resuspend the beads.
9. Place the deepwell plate on the magnetic support and remove the supernatant without touching the magnetic beads.
10. Remove the plate from the support, add 200  $\mu$ L of WASH 3 and resuspend the beads.
11. Place the deepwell plate on the magnetic support and remove the supernatant and any drops without touching the magnetic beads.
12. Remove the plate from the support, add 120  $\mu$ L of Elution Buffer and resuspend the beads.
13. Incubate for 5 minutes at room temperature.
14. Place the deepwell plate on the magnetic support and collect the eluate in a new plate or tube.

**Note:** For short-term storage of less than 24 hours, we recommend storing the eluate at +4°C. For longer storage, the eluate should be stored at -80°C.

## Quality control

The products in the kit are tested using the above protocol on an artificial matrix. The extraction product is then quantified and compared with our reference data.