Instructions d'utilisation - Kit 96-24-14- 1réaction(s) (p.1-2)**Instructions for use – 96-24-14-1 reaction(s) kit (p.3-4)**

Nucleis Mag Extractor LIBERTY

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Nucleis Mag Extractor LIBERTY

Shipping and storage conditions

Nucleis Mag Extractor LIBERTY is shipped at room temperature. On arrival, the kit should be stored at +4°C. These components remain stable until the indicated expiry date.

Composition

This kit, for complex matrices such as serum, plasma, urine, synovial fluid, tissue, stool, amniotic fluid, bone marrow, bronchoalveolar lavage (BAL), nasopharyngeal swab, ocular swab, urethral swab, and cutaneous-mucosal sample., enables up to 96 samples to be extracted from an initial sample volume of 150 µL.

The different products contained in the kit are described in the table below:

Components	Volumes for 96 reactions
PK (Proteinase K)	1.1 mL
RNA Carrier	0.12 mL
PVG	18 mL
SP	6 mL
BMB	15 mL
WASH 1	22 mL
WASH 2	22 mL
WASH 3	22 mL
Tampon d'élution	15 mL

Components	Volumes for 24 reactions	Volumes for 14 reactions	Volumes for 1 reactions
PK (Proteinase K)	0.3 mL	0.2 mL	10 µL
PVG + RNACarrier	4.5 mL	2.4 mL	180 µL
SP	1.5 mL	0.8 mL	100 µL
BMB	4 mL	2 mL	150 µL
WASH 1	5 mL	3.5 mL	225 µL
WASH 2	5 mL	3.5 mL	225 µL
WASH 3	5 mL	3.5 mL	225 µL
Tampon d'élution	4 mL	2 mL	180 µL

Equipment required

For manual use of this kit, the equipment (not supplied in the kit) listed below is required:

- Adjustable 10 µL and 200 µL micropipettes
- Sterile low retention filter cones for micropipettes
- Deepwell 96 wells
- Adhesive aluminium foil
- Incubator up to 56°C
- Magnetic support for deepwell 96 wells

Manual protocol

Depending on the matrix, it is sometimes necessary to carry out a pre-treatment before starting the extraction. The file listing the pre-treatments for each matrix can be found on our website.

1. In a deepwell plate, insert 150 μ L of sample, add 150 μ L of PVG + 1 μ L RNACarrier and 6 μ L of PK. Mix by pipetting gently.
2. Cover the plate with adhesive aluminium foil and incubate for 35 minutes at 42°C.
3. Remove the aluminium foil. Add 40 μ L of SP and 125 μ L of BMB, mix by pipetting gently.
4. Incubate for 5 minutes at room temperature.
5. Place the deepwell plate on the magnetic support and remove the supernatant without touching the magnetic beads.
6. Remove the plate from the support, add 200 μ L of WASH 1 and resuspend the beads.
7. Place the deepwell plate on the magnetic support and remove the supernatant without touching the magnetic beads.
8. Remove the plate from the support, add 200 μ L of WASH 2 and resuspend the beads.
9. Place the deepwell plate on the magnetic support and remove the supernatant without touching the magnetic beads.
10. Remove the plate from the support, add 200 μ L of WASH 3 and resuspend the beads.
11. Place the deepwell plate on the magnetic support and remove the supernatant and any drops without touching the magnetic beads.
12. Remove the plate from the support, add 120 μ L of Elution Buffer and resuspend the beads.
13. Incubate for 5 minutes at room temperature.
14. Place the deepwell plate on the magnetic support and collect the eluate in a new plate or tube.

Note: For short-term storage of less than 24 hours, we recommend storing the eluate at +4°C. For longer storage, the eluate should be stored at -80°C.

Quality control

The products in the kit are tested using the above protocol on an artificial matrix. The extraction product is then quantified and compared with our reference data.